

**The Nottingham Total Laparoscopic Hysterectomy  
Masterclass Course  
Utilising Cadaveric and Clinical Immersion  
Training Model**

**Course Director: Mr Jafaru Abu**

**Cadaver Lab Learning Targets:**

1. Understand how to use the following instruments:
  - i. Ligasure
  - ii. McCartney tube
2. Safe pneumoperitoneum
3. Place ports appropriately for good ergonomics
4. Able to perform Palmer's pneumoperitoneum and initial entry
5. Technique of Hasson's entry
6. Expose pelvic sidewall
7. Expose ureter, external iliac vessels, pararectal/paravesical spaces
8. Secure IP ligament with direct ureteric visualisation
9. Safely reflect the bladder away from cervix and upper vagina
10. Safely reflect the rectum away from cervix and upper vagina
11. Perform posterior colpotomy, extend laterally on both sides
12. Perform anterior colpotomy
13. Close the vaginal vault with the v-loc
14. Close the vaginal using O-monocryl
15. Make a hole in the bladder and close with 2/0 monocryl tying knots laparoscopically